Report No.（报告编号）：20260826-912-01

Test report

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Experimental Requirements 实验要求
Separation between the main substance and custom peptide shall be greater than 1.5; signal-to-noise ratio for sensitivity test (limit of 0.05%) shall be greater than 10; peak tailing shall be improved. 要求主物质与定制肽之间的分离度＞1.5，灵敏度（0.05%限度）信噪比＞10，改善拖尾

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1.2. Solution Preparation（溶液配制）

1.2.1. Mobile Phase Preparation（流动相配制）：

Pipette 0.5 mL trifluoroacetic acid, dilute with ultrapure water to 500 mL, then add 500 mL acetonitrile, mix well and degas by sonication.

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1.2.3. Mixed solution（混合溶液）：

Weigh 0.5 mg each of the customer-provided main substance and custom peptide, dissolve with mobile phase to prepare a mixed solution containing 0.5 mg/mL of each component, and shake well.

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Weigh 1 mg of the main substance supplied by the customer, dissolve with mobile phase to prepare a solution of 1 mg/mL, shake well to obtain the solution.

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Weigh 1 mg of the custom peptide provided by the customer, dissolve with mobile phase to prepare a 1 mg/mL solution, and shake well to obtain the single standard solution.称取客户提供的定制肽 1mg，加流动相溶解制成 1mg/ml 的溶液，摇匀，即得。

1.2.7. Sensitivity solution（灵敏度溶液）：

Take an appropriate volume of single standard solution of the main substance, dilute with mobile phase to prepare a solution at 0.5 μg/mL, shake well to obtain the solution.取主物质单标溶液适量，加流动相稀释制成 0.5 μg/ml 的溶液，摇匀，即得。

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